



US FDA CDER DMF No.: 038447 US FDA CBER DMF No.: 29657 Acute Systemic Toxicity Test

INTRODUCTION

DASEA® Ultramedia® CIK cell serum-free medium is speially designed to promote high-efficiency proliferation of CIK cells in vitro. DASEA® Ultramedia® CIK contains no animal-derived components or heterologous growth factors which enhances experimental reproducibility by avoiding batch-to-batch variations caused by serum.

* For Research Use Only



High Proliferation Efficiency



High Positive Rate



GMP Manufacturing

ADVANTAGES

Safe and reliable raw material management



- Thousand-square-meter storage space for raw materials.
- Strictly monitor temperature and humidity of the warehouse.



Optimal batch-to-batch stability

• Compliance with raw material release criteria.

マ 中科睿极

产品货号: RGL0040

STERILE A

DASEA® Ultramedia® CIK 细胞无血清培养基

2°C 18°C

本产品仅用于科研用途


• Strictly define and enforce production process control parameters.



QUALITY STANDARD


Quality Inspection Item	Quality Inspection Standards
Appearance	Orange-red clear liquid
рН	7.2 ± 0.2
Osmolality	280-320 mOsm/kg
Endotoxin	<0.5 EU/mL
Sterility Test	Negative
Mycoplasma Test	Negative

PERFORMANCE DATA



D0: Total 10ml blood was used to collect CIK cells. Initial cell culture volume is 5mL and initial cell count is 7.5×10^6 .

D14: Cell number were increased by 500–1000 times. 3.7 billion to 7.5 billion cells were harvested, with 1.2–2L DASEA Ultramedia $^{\circ}$ CIK medium usage.



The cell viability kept high during the high density cell culture process.

