

DASEA® UltraTryple® recombinant enzymes

FDA CDER DMF: 039384 FDA CBER DMF: 30280



PRODUCT INTRODUCTION

UltraTryple® is a recombinant serine protease free of animal-derived components. Compared to traditional trypsin, it has an extremely low risk of contamination from exogenous viruses, offering high biosafety. It can completely replace porcine or bovine trypsin in various fields of bioengineering. It is commonly used for applications such as adherent cell digestion, primary cell isolation, and the production of recombinant proteins and vaccines. Utilizing a calcium- and magnesium-free PBS buffer system, it provides mild digestion suitable for cell passaging processes in both serum-free and serum-containing systems. This product contains DPBS/1.1mM EDTA, which are stable at room temperature, and easy to store and use.


ADVANTAGES


- **Gentle**
 - Ensures minimal cell damage, even after 30 minutes of consecutive digestion, with no effect on the phenotype and proliferation of protease-sensitive cells;
- High Efficiency and Convenience
 No inhibitor required, easy termination by dilution;
- Extensive Application
 Suitable for a diverse range of cell types (MSC, 293, VERO, MDCK, etc.);
- Time-Saving
 No thawing required, reduced preparation time, rapid cell separation, short digestion time;
- Animal Origin-Free (AOF)
 No animal-derived components;
 - Stability
 Long-term storage at room temperature, no refrigeration or -20°C storage required for enhanced convenience.

PERFORMANCE DATA




UltraTryple® outperforms Competitor 1 (recombinant) and Competitor 2 (animal-derived) in terms of fold expansion after consecutive subculturing of stem cells, with each passage's proliferation efficiency increasing by 5-20%.








^{*} Recover P3 of MSCs from cryogenic storage, Seed them at a density of 8,000 cells per square centimeter (cm²) in T75 flasks. After consecutive subculture for 3 passages, use UltraTryple*, Competitor 1 (recombinant) and Competitor 2 (animal-derived) to digest the cells respectively.

Under both normal and extended (30min) digestion conditions, UltraTryple® outperforms Competitor 1 (recombinant) and Competitor 2 (animal-derived) in preserving cell morphology.









^{*} After three consecutive subcultures and 3 minutes of digestion, Competitor 2 digestion resulted in a decreased percentage of CD105-positive cells, while UltraTryple* digestion maintained normal cell morphology.

^{*} After digestion for more than 30 minutes, Competitor 2 digestion resulted in a significantly decreased percentage of CD105-positive cells, while UltraTryple® digestion maintained normal cell morphology.



UltraTryple® is suitable for a diverse range of cell types, including strongly adherent cells. Moreover, it offers faster digestion times compared to competitors.

| | Competitor 1 (recombinant) | | Competitor 2 (animal-derived) | | UltraTryple® | |
|-----------|-------------------------------|-----------------------|----------------------------------|-----------------------|---------------------|-----------------------|
| Cell Type | ls it digestible | Digestive duration | ls it digestible | Digestive duration | ls it digestible | Digestive duration |
| MSC | Easy | 3m30s | Easy | 3m30s | Easy | 3min |
| 293 | Easy | 7m30s | Easy | Twice | Easy | 6m30s |
| MDCK | Easy | 50min | Difficult | / | Easy | 40min |
| VERO | Easy | 15min | Difficult | / | Easy | 13min |

^{*} Test data is for reference only

| Item | Competitor 1 (recombinant) | Competitor 2 (animal-derived) | Ultra Tryple® |
|-----------------|-------------------------------|---|--------------------|
| Storage | Normal | −20°C | Normal |
| Transportation | Normal | Wet ice | Normal |
| Termination | Supernatant or PBS | Serum, Complete culture medium or protease inhibitor | Supernatant or PBS |
| Virus Detection | Not required | Required | Not required |
| Animal-derived | None | Contain | None |
| Effect | None | Affects surface proteins | None |
| EDTA | 1.1mM | 0.9mM | 1.1mM |
| Capacity | **** | ★★★ ☆☆ | **** |
| Concentration | 1x | 0.25%, can be diluted to 0.05% | Customizable |

QUALITY STANDARD

| Inspection Items | Inspection Standards | | |
|-------------------|----------------------|--|--|
| Sterility Test | Not detected | | |
| Osmolality | 260-340 mOsm/kg | | |
| рН | 7.0 - 7.4 | | |
| Endotoxin | <0.25EU/mL | | |
| In Vitro Bioassay | Qualified | | |
| Mycoplasma | Not detected | | |
| Appearance | Pale yellow liquid | | |

Enterprise Quality Management System









ORDER INFORMATION

| Product | Package (per bottle) | P/N | Storage | Shelf Life |
|---------------------------|----------------------|---------|-----------------------|------------|
| DASEA® UltraTryple® | 100mL | RGM0061 | Protect from light | 1 year |
| recombinant enzymes | 500mL | RGM0062 | | |
| DASEA® UltraTryple® | 100mL | RGM0064 | | |
| recombinant enzymes (10X) | 500mL | RGM0065 | | |

Related Products

| | Recombinant enzymes 100mL*2 units PBS 500mL*2 | RGE0071 | | |
|---|---|---------|-----------------------|--------|
| Recombinant enzymes kit | Recombinant enzymes (10X) 100mL*2 units PBS 500mL*2 units EDTA (0.5M) 1mL*2 units | RGE0072 | Protect from light | 1 year |
| | Recombinant enzymes 500mL*10 units PBS 500mL*10 units | RGE0073 | | |
| PBS | 500mL/bottle | RGE0031 | | |



REGEN- α GEEK (Haining) Biotechnology Co., Ltd.

Tel: +860755-26412015

Web: http://www.regengeek.com Email: enquiry@regengeek.com
Address: No. 301, Building 6, Juanhu Science and Innovation Park, No. 500, Shuiyue Pavilion East
Road, Xiashi Street, Haining City, Jiaxing City, Zhejiang Province, China.