マ。 中科睿极

DASEA* Ultramedia* NK细胞无血清培养基



2'0 1 8'0 DASEA® Ultramedia® NK cell induction and expansion culture kit

US FDA CDER DMF No.: 038393 **US FDA CBER DMF No.: 29527 Acute Systemic Toxicity Test**

INTRODUCTION

DASEA® Uitramedia® NK cell induction and expansion culture kit is a serum-free culture medium developed by REGEN-αGEEK, that contains no heterologous components and has a clear composition. It is suitable for stimulating activation and continuous expansion of NK cells in Peripheral Blood Mononuclear Cells (PBMCs) isolated from human peripheral blood.

This kit uses a pure factor method and does not contain feeder cells, avoiding ethical issues and residual oncogenic risks associated with feeder cells. The components of this product are clear and quality-controlled, and its production follows the ISO9001 quality system and GMP guidelines.

* For research use only

ADVANTAGES



High proliferation rate



No heterologous Components



Pure factor method



The median of cell positivity rate is greater than 70%




Easy to operate



manufacturing

MORPHOLOGY







During the culture process, the number of cell clusters increased, and cell morphology is uniform and the culture appears transparent and bright.

PRODUCT COMPOSITION

Product Number	Description	Package	Shelf Life	Storage Conditions
RGL0011	NK cell serum-free medium	1L	12 months	Protect from light / 2°C~8°C
	NK factor A	50 μL		-25°C ~ -10°C
	NK factor B	1 unit		-25°C ~ -10°C

PERFORMANCE DATA







Compared with competitors, the NK cell induction and expansion culture kit can amplify NK cells from PBMCs more efficienntly.

In 1L cell culture system, the cell expansion fold was 150–240 times, and the number of harvested cells was $20x10^8$ ~36x10⁸.

Note: The actual number of harvested cells is closely related to the quality of cells themselves.



The positive rate of samples in the figure is the summary of test data, and the median positive rate of cells is more than 70%.



The cytotoxic efficiency of NK against K562 at Day 14 under different effect target ratios.

Cytotoxicity kit: CFSE/7-AAD

